

2'-O-methylcytidine/ Cm Rabbit mAb

Catalog No.: A20693 **Recombinant**

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

DB, ELISA

Cross-Reactivity

Species independent

CloneNo number

ARC50639

Background

RNA methylation plays a significant regulatory role in various of physiological activities. 2'-O-methylcytidine (Cm) is one of the representative 2'-O-methylation and base-methylation modified epigenetic marks of RNA. 2'-O-Methylcytidine is derivative of the nucleoside residue cytidine with a methyl group attached to the 2'-oxygen position. 2'-O-methylcytidine was found to be related to various diseases including cancer, and it might have great potential to be novel biomarkers for detection of breast cancer in the early stage. In addition, 2'-O-methylcytidine showed as an inhibitor of RNA polymerase from the hepatitis C virus (HCV). And 2'-O-methylcytidine was involved in plant stress responses.

Recommended Dilutions

DB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Chemical compounds corresponding to 2'-O-methylcytidine/ Cm.

Synonyms

Cm; 2'-O-methylcytidine; 2'-O-methylcytidine/ Cm

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol and 0.05% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



The 2'-O-methylcytidine(Cm) Rabbit mAb (A20693) are tested in Dot Blot against 2'-O-methylcytidine and unmodified cytidine.2'-O-methylcytidine
□Biotin-5'CGATAACCACTAGT(Cm)3'unmodified cytosine□Biotin-5'CGATAACCACTAGTC3'