

SLIT3 Rabbit pAb

Catalog No.: A20741

Basic Information

Observed MW

Refer to figures

Calculated MW

168kDa

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

The protein encoded by this gene is secreted, likely interacting with roundabout homolog receptors to effect cell migration. Two transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

WB 1:500 - 1:1000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

6586

Swiss Prot

O75094

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

MEGF5; SLIL2; SLIT1; slit2; Slit-3; SLIT3

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.