

USP34 Rabbit pAb

Catalog No.: A20744

Basic Information

Observed MW

Refer to figures

Calculated MW

404kDa

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Enables cysteine-type endopeptidase activity and thiol-dependent deubiquitinase. Involved in positive regulation of canonical Wnt signaling pathway and protein K48-linked deubiquitination. Predicted to be active in cytosol and nucleus.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

9736

Swiss Prot

Q70CQ2

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.