

ATFS-1Rabbit pAb

Catalog No.: A20751

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Caenorhabditis elegans

Background

Enables DNA-binding transcription factor activity, RNA polymerase II-specific; enzyme binding activity; and mitochondrial transcription factor activity. Involved in several processes, including mitochondrial unfolded protein response; positive regulation of nematode larval development; and regulation of gene expression. Located in mitochondrial matrix and nucleus. Is expressed in intestine.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

179922

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.