

PSTPIP2 Rabbit pAb

Catalog No.: A20765

Basic Information

Observed MW

Refer to figures

Calculated MW

39kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable actin filament binding activity. Predicted to be involved in actin filament polymerization. Predicted to be located in cytoskeleton and membrane. Predicted to be active in actin filament; cytoplasm; and plasma membrane.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

9050

Swiss Prot

Q9H939

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

MAYP; PSTPIP2

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.