

ADGRG3 Rabbit pAb

Catalog No.: A20780

Basic Information

Observed MW

Refer to figures

Calculated MW

61kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable G protein-coupled receptor activity. Involved in G protein-coupled receptor signaling pathway and regulation of cell migration. Located in plasma membrane.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

222487

Swiss Prot

Q86Y34

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

PB99; GPR97; PGR26; ADGRG3

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.