

RABGEF1 Rabbit pAb

Catalog No.: A20789 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

68kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

RABGEF1 forms a complex with rabaptin-5 (RABPT5; MIM 603616) that is required for endocytic membrane fusion, and it serves as a specific guanine nucleotide exchange factor (GEF) for RAB5 (RAB5A; MIM 179512) (Horiuchi et al., 1997 [PubMed 9323142]).

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

27342

Swiss Prot

Q9HCQ5

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 222-491 of human RABGEF1 (NP_055319.1).

Synonyms

RAP1; RABEX5; rabex-5; RABGEF1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.