

BrdU Rabbit mAb

Catalog No.: A20790

Recombinant

2 Publications

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

ELISA,IF/ICC

Cross-Reactivity

Species independent

CloneNo number

ARC51030

Background

Recommended Dilutions

IF/ICC 1:50 - 1:200

Immunogen Information

Gene ID

Swiss Prot

Immunogen

BrdU

Synonyms

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.

Validation Data



Immunofluorescence analysis of BrdU treated (left figure) and untreated (right figure) U2OS cells using BrdU antibody (A20790) at dilution of 1:100. Blue: DAPI for nuclear staining.