

TOE1 Rabbit pAb

Catalog No.: A20801

Basic Information

Observed MW

Refer to figures

Calculated MW

57kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Enables poly(A)-specific ribonuclease activity and snRNA binding activity. Involved in RNA phosphodiester bond hydrolysis, exonucleolytic and snRNA 3'-end processing. Located in Cajal body and cytoplasm. Implicated in pontocerebellar hypoplasia type 7.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

114034

Swiss Prot

Q96GM8

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

PCH7; TOE-1; hCaf1z; TOE1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.