

HB1 Rabbit pAb

Catalog No.: A21284

Basic Information

Observed MW

Refer to figures

Calculated MW

31kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Arabidopsis thaliana

Background

Encodes a homeodomain leucine zipper class I (HD-Zip I) transcriptional activator involved in leaf development.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

821138

Swiss Prot

Q02283

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ARABIDOPSIS THALIANA HOMEODOMAIN PROTEIN 1; ATHB-1; ATHB1; F4P13.2; F4P13_2; HAT5; HD-ZIP-1; homeobox 1; HOMEODOMAIN PROTEIN FROM ARABIDOPSIS THALIANA 5; HB1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.