

TTC3 Rabbit pAb

Catalog No.: A21639

Basic Information

Observed MW

Refer to figures

Calculated MW

230kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables ubiquitin-protein transferase activity. Involved in protein K48-linked ubiquitination and ubiquitin-dependent protein catabolic process. Located in cytosol; nucleolus; and nucleoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

7267

Swiss Prot

P53804

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

DCRR1; RNF105; TPRDIII; TTC3

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.