

VRN1 Rabbit pAb

Catalog No.: A22256

Basic Information

Observed MW

Calculated MW

39kDa

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Arabidopsis thaliana

Background

Required for vernalization. Essential for the complete repression of FLC in vernalized plants. Required for the methylation of histone H3

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

821432

Swiss Prot

Q8L3W1

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

REDUCED VERNALIZATION RESPONSE 1; REM39; REPRODUCTIVE MERISTEM 39; VRN1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.