

# CASS4 Rabbit pAb

Catalog No.: A22559

## Basic Information

### Observed MW

### Calculated MW

87kDa

### Category

Primary antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human, Mouse, Rat

## Background

Enables protein tyrosine kinase binding activity. Involved in several processes, including positive regulation of protein kinase B signaling; positive regulation of protein tyrosine kinase activity; and positive regulation of substrate adhesion-dependent cell spreading. Located in focal adhesion. Part of cytoplasm.

## Recommended Dilutions

**WB** 1:500 - 1:2000

**ELISA** Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

57091

### Swiss Prot

Q9NQ75

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

CAS4; HEFL; HEPL; C20orf32; CASS4

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.