

Rabbit Anti-Mouse IgG2a

Catalog No.: A22752

Basic Information

Observed MW

50-55kDa

Calculated MW

36kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Mouse

CloneNo number

ARC57519

Background

Recommended Dilutions

WB 1:10000 - 1:190000

Immunogen Information

Gene ID

□

Swiss Prot

P01863

Immunogen

Recombinant protein of human Rabbit Anti-Mouse IgG2a.

Synonyms

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.

Validation Data



Western blot analysis of various lysates, using Rabbit Anti-Mouse IgG2a (A22752) at 1:160000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25 μ g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.