

ZmPUB32 Rabbit pAb

Catalog No.: A22856

Basic Information

Observed MW

Calculated MW

90kDa

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Maize

Background

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

100382053

Swiss Prot

C0P3Y9

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

GRMZM2G101754; ZmPUB32

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.