

ATG16L2 Rabbit pAb

Catalog No.: A22862

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse

Background

Acts upstream of or within negative stranded viral RNA replication. Located in cytosol. Is expressed in embryo. Orthologous to human ATG16L2 (autophagy related 16 like 2).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

73683

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

2410118P20Rik; ATG16L2

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.