

# FITC Mouse IgG2a isotype control

Catalog No.: A22880

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

### Category

Primary antibody

### Applications

FC

### Cross-Reactivity

### CloneNo number

AMC0554

### Conjugate

FITC. Ex:491nm. Em:516nm.

## Background

---

## Recommended Dilutions

---

FC 10 µL/Test

## Immunogen Information

---

Gene ID

Swiss Prot

Immunogen

Synonyms

## Contact

---

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

---

### Source

Mouse

### Isotype

Mouse IgG2a

### Purification

Affinity purification


### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.1% Proclin300,0.2% BSA,pH7.2.

## Validation Data

---



Flow cytometry:  $1 \times 10^6$  Human PBMC were surface-stained with FITC Mouse anti-Human CD235a mAb (A22831, orange line) or FITC Mouse IgG2a isotype control mAb (A22880, blue line).