

APC Mouse IgG2b isotype control

Catalog No.: A22884

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

FC

Cross-Reactivity

CloneNo number

AMC0555

Conjugate

APC. Ex:650nm. Em:660nm.

Background

Recommended Dilutions

FC 10 μ L/Test

Immunogen Information


Gene ID

Swiss Prot

Immunogen

Synonyms

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Mouse

Isotype

Mouse IgG2b

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.1% Proclin300,0.2% BSA,pH7.2.

Validation Data



Flow cytometry: 1×10^6 Human PBMC were surface-stained with APC Mouse anti-Human CD11b mAb (A22811, orange line) or APC Mouse IgG2b isotype control mAb (A22884, blue line).