

TPPP3 Rabbit mAb

Catalog No.: A22991 **Recombinant**

Basic Information

Observed MW

19kDa

Calculated MW

19kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Mouse, Rat

CloneNo number

ARC58837

Background

Enables tubulin binding activity. Involved in decidualization and microtubule bundle formation. Colocalizes with microtubule bundle and perinuclear region of cytoplasm.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

51673

Swiss Prot

Q9BW30

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-176 of human TPPP3 (NP_057048.2).

Synonyms

p20; CGI-38; TPPP/p20; p25gamma; TPPP3

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using TPPP3 Rabbit mAb (A22991) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25 µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 60s.