

Goat IgG (Fc) fragment Rabbit mAb

Catalog No.: A23074 **Recombinant**

Basic Information

Observed MW

62kDa

Calculated MW

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Goat

CloneNo number

ARC56943

Background

Recommended Dilutions

WB 1:2000 - 1:6000

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Goat IgG

Synonyms

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300, 0.05% BSA, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of recombinant Goat IgG (Fc) fragment and Goat IgG (Fab) fragment [Negative control], using Goat IgG (Fc) fragment Rabbit mAb (A23074) at 1:5000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 10ng per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 180s.