

# PE Anti-Mouse CD19 mAb

Catalog No.: A23331

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

### Category

Primary antibody

### Applications

FC

### Cross-Reactivity

Mouse

### CloneNo number

ATC0004

### Conjugate

PE. Ex:565nm. Em:574nm.

## Recommended Dilutions

**FC** The amount of antibody required for optimal staining of a cell sample should be determined empirically in your system.

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Background

Involved in several processes, including B-1 B cell differentiation; positive regulation of phosphatidylinositol 3-kinase activity; and positive regulation of release of sequestered calcium ion into cytosol. Acts upstream of or within B cell receptor signaling pathway. Located in external side of plasma membrane. Is integral component of plasma membrane. Is expressed in liver and spleen. Human ortholog(s) of this gene implicated in common variable immunodeficiency. Orthologous to human CD19 (CD19 molecule).

## Immunogen Information

### Gene ID

12478

### Swiss Prot

### Immunogen

Recombinant protein of Mouse CD19.

### Synonyms

## Product Information

### Source

Rat

### Isotype

Rat IgG2a, kappa

### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: 10 mM NaH<sub>2</sub>PO<sub>4</sub>, 150 mM NaCl, 0.09% NaN<sub>3</sub>, 0.1% gelatin, pH7.2

## Validation Data

---

