

DNAJC4 Rabbit pAb

Catalog No.: A23874

Basic Information

Observed MW

Calculated MW

28kDa

Category

Primary antibody

Applications

ELISA,IF/ICC

Cross-Reactivity

Human, Rat

Background

Predicted to be involved in response to unfolded protein. Predicted to be integral component of membrane.

Recommended Dilutions

IF/ICC 1:50 - 1:200

Immunogen Information

Gene ID

3338

Swiss Prot

Q9NNZ3

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 35-100 of human DNAJC4(NP_005519.2).

Synonyms

HSPF2; MCG18; DANJC4; DNAJC4

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.

Validation Data



Immunofluorescence analysis of HepG2 cells using DNAJC4 Rabbit pAb (A23874) at dilution of 1:20 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of PC-12 cells using DNAJC4 Rabbit pAb (A23874) at dilution of 1:20 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.