

ABflo® 647 Rabbit anti-Human/Mouse CD55/DAF mAb

Catalog No.: A23898

Basic Information

Observed MW

Calculated MW

42KDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human, Mouse

CloneNo number

ARC60934-ABflo647

Conjugate

ABflo® 647. Ex:648nm. Em:664nm.

Recommended Dilutions

FC 5 µl per 10⁶ cells in
100 µl volume

Background

This gene encodes an inhibitor of both the classical and the alternative pathways of complement activation. The encoded preproprotein undergoes post-translational processing to generate a mature polypeptide anchored to the plasma membrane via a glycosylphosphatidylinositol moiety. Erythrocytes from mice deficient in the encoded protein exhibit impaired regulation of complement activation resulting in enhanced complement deposition. Mice lacking the encoded protein exhibit enhanced susceptibility to experimentally induced myasthenia gravis. This gene is located adjacent to a closely related gene on chromosome 1.

Immunogen Information

Gene ID

13136

Swiss Prot

Q61475

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 35-361 of mouse CD55/DAF [NP_034146.2].

Synonyms

Daf; Daf1; Daf-GPI; GPI-DAF

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.03% proclin300, 0.2% BSA, pH7.3.


Validation Data




Flow cytometry: 1×10^6 C57BL/6 mouse Splenocytes were surface-stained with ABflo® 647 Rabbit anti-Human/Mouse CD55/DAF mAb(A23898,5 μ l/Test,orange line) or ABflo® 647 Rabbit IgG isotype control (A22070,5 μ l/Test,blue line).Non-fluorescently stained C57BL/6 mouse Splenocytes were used as blank control (red line).



Flow cytometry: 1×10^6 SH-SY5Y cells(negative control,left) and Hela cells(right) were surface-stained with ABflo® 647 Rabbit anti-Human/Mouse CD55/DAF mAb(A23898,5 μ l/Test,orange line) or ABflo® 647 Rabbit IgG isotype control (A22070,5 μ l/Test,blue line).Non-fluorescently stained cells were used as blank control (red line). ; Flow cytometry: 1×10^6 C57BL/6 mouse Splenocytes were surface-stained with ABflo® 647 Rabbit anti-Human/Mouse CD55/DAF mAb(A23898,5 μ l/Test,orange line) or ABflo® 647 Rabbit IgG isotype control (A22070,5 μ l/Test,blue line).Non-fluorescently stained C57BL/6 mouse Splenocytes were used as blank control (red line).



Flow cytometry: 1×10^6 C57BL/6 mouse Splenocytes were surface-stained with ABflo® 647 Rabbit IgG isotype control (A22070,5 μ l/Test,left) or ABflo® 647 Rabbit anti-Human/Mouse CD55/DAF mAb(A23898,5 μ l/Test,right).



Flow cytometry: 1×10^6 Hela cells were surface-stained with ABflo® 647 Rabbit IgG isotype control (A22070,5 μ l/Test,left) or ABflo® 647 Rabbit anti-Human/Mouse CD55/DAF mAb(A23898,5 μ l/Test,right).