

RINL Rabbit pAb

Catalog No.: A23943

Basic Information

Observed MW

Refer to figures

Calculated MW

62kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables guanyl-nucleotide exchange factor activity. Acts upstream of or within endocytosis. Located in actin cytoskeleton and ruffle. Orthologous to human RINL (Ras and Rab interactor like).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

320435

Swiss Prot

Q80UW3

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

5830482F20Rik; 9930116N10Rik; RINL

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.