Smad9 Rabbit pAb

Catalog No.: A24036



Basic Information

Observed MW

Calculated MW

52kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

The protein encoded by this gene is a member of the SMAD family, which transduces signals from TGF-beta family members. The encoded protein is activated by bone morphogenetic proteins and interacts with SMAD4. Two transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

Immunogen Information

Gene ID Swiss Prot4093
015198

Immunogen

 $Synthetic\ peptide.\ This\ information\ is\ considered\ to\ be\ commercially\ sensitive.$

Synonyms

PPH2; MADH6; MADH9; SMAD8; SMAD8A; SMAD8B; SMAD8/9; Smad9

Contact

	400-999-6126
×	cn.market@abclonal.com.cn
	www.abclonal.com.cn

Product Information

SourceIsotypePurificationRabbitIgGAffinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.