

# AQP4 Rat pAb

Catalog No.: A24050

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

34kDa

### Category

Primary antibody

### Applications

WB,IHC-P,IF/IHC-Fr,ICC

### Cross-Reactivity

Mouse, Rat, Human

## Background

## Recommended Dilutions

<b>WB</b>	1:1000-2500
<b>IHC-P</b>	1:500-1000
<b>ICC</b>	1:1000
<b>IF/IHC-Fr</b>	1:500-2000

## Immunogen Information

### Gene ID

361

### Swiss Prot

P55087

### Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

### Synonyms

AQP4

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rat

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze/thaw cycle

Buffer:pH7.4, PBS, 55mM Tris, 50mM Glycine, 40% Glycerol, 0.1% potassium sorbate