

MYRF Rabbit pAb

Catalog No.: A24051

Basic Information

Observed MW

Refer to figures

Calculated MW

75kDa

Category

Primary antibody

Applications

IHC-P,IF/IHC-Fr,ICC

Cross-Reactivity

Mouse, Rat, Human, Chicken, Dog

Background

Recommended Dilutions

IHC-P	1:200-500
ICC	1:500
IF/IHC-Fr	1:500-1000

Immunogen Information

Gene ID

225908

Swiss Prot

Q3UR85

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

MYRF

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze/thaw cycle

Buffer: pH7.4, PBS, 55mM Tris, 50mM Glycine, 40% Glycerol, 0.1% potassium sorbate