

Ctip2 Guinea pig pAb

Catalog No.: A24052

Basic Information

Observed MW

Refer to figures

Calculated MW

95kDa

Category

Primary antibody

Applications

WB,IHC-P,IF/IHC-Fr,ICC

Cross-Reactivity

Mouse, rat, human

Background

Recommended Dilutions

WB	1:2000-5000
IHC-P	1:500-1000
ICC	1:1000
IF/IHC-Fr	1:1000-2000

Immunogen Information

Gene ID

58208

Swiss Prot

Q99PV8

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

Ctip2

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Guinea pig

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C.Avoid freeze/thaw cycle

Buffer:pH7.4, PBS, 55mM Tris, 50mM Glycine, 40% Glycerol, 0.1% potassium sorbate