

CNPase Rat pAb

Catalog No.: A24053

Basic Information

Observed MW

Refer to figures

Calculated MW

47kDa

Category

Primary antibody

Applications

WB,IHC-P,IF/IHC-Fr,ICC

Cross-Reactivity

Mouse, Rat, Human

Background

Recommended Dilutions

WB	1:1000-2000
IHC-P	1:100-250
ICC	1:250
IF/IHC-Fr	1:100-250

Immunogen Information

Gene ID

1267

Swiss Prot

P09543

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

CNPase

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rat

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze/thaw cycle

Buffer:pH7.4, PBS, 55mM Tris, 50mM Glycine, 40% Glycerol, 0.1% potassium sorbate