

Cytokeratin(HMW) Mouse mAb

Catalog No.: A24141

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

IHC-P

Cross-Reactivity

Human

CloneNo number

AMC0667

Background

Recommended Dilutions

IHC-P 1:100 - 1:200

Immunogen Information

Gene ID

Swiss Prot

Immunogen

This information is considered to be commercially sensitive.

Synonyms

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Mouse

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol and 0.05% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Immunohistochemistry analysis of paraffin-embedded Human prostate tissue using Cytokeratin(HMW) Mouse mAb (A24141) at a dilution of 1:100 (20x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.