ABclonal® www.abclonal.com

Neuron Specific Enolase(NSE) Mouse mAb

Catalog No.: A24147

Basic Information

Observed MW

Calculated MW

47kDa

Category

Primary antibody

Applications

IHC-P

Cross-Reactivity

Human

CloneNo number

AMC0673

Background

This gene encodes one of the three enolase isoenzymes found in mammals. This isoenzyme, a homodimer, is found in mature neurons and cells of neuronal origin. A switch from alpha enolase to gamma enolase occurs in neural tissue during development in rats and primates.

Recommended Dilutions

IHC-P

1:100 - 1:200

Immunogen Information

Gene ID 2026 **Swiss Prot**

P09104

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

NSE; HEL-S-279

Contact

a	400-999-6126
×	cn.market@abclonal.com.cn
	www.abclonal.com.cn

Product Information


SourceIsotypePurificationMouseIgGAffinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol and 0.05% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Immunohistochemistry analysis of paraffinembedded Human renal cell carcinoma tissue using Neuron Specific Enolase(NSE) Mouse mAb (A24147) at a dilution of 1:100 (20x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.