

ABflo® 594 Rabbit anti-Human CD72 mAb

Catalog No.: A24238

Basic Information

Observed MW

Calculated MW

40kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human

CloneNo number

ARC62240-ABflo594

Conjugate

ABflo® 594. Ex:588nm. Em:604nm.

Background

Predicted to enable signaling receptor binding activity. Predicted to be involved in cell adhesion. Is integral component of plasma membrane.

Recommended Dilutions

FC 5 µl per 10⁶ cells in
100 µl volume

Immunogen Information

Gene ID

971

Swiss Prot

P21854

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 117-359 of human CD72(NP_001773.1).

Synonyms

CD72; CD72b; LYB2; CD72 molecule

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.03% proclin300,0.2% BSA,pH7.3.

Validation Data



Flow cytometry: 1×10^6 A549 cells (negative control, left) and Daudi cells (right) were surface-stained with ABflo® 594 Rabbit anti-Human CD72 mAb (A24238, 5 μ l/Test, orange line) or ABflo® 594 Rabbit IgG isotype control (A23821, 5 μ l/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).



Flow cytometry: 1×10^6 Daudi cells were surface-stained with ABflo® 594 Rabbit IgG isotype control (A23821, 5 μ l/Test, left) or ABflo® 594 Rabbit anti-Human CD72 mAb (A24238, 5 μ l/Test, right).