

ABflo® 594 Rabbit anti-Human CD112/Nectin-2 mAb

Catalog No.: A24329

Basic Information

Observed MW

Refer to figures

Calculated MW

51kDa/57kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human

CloneNo number

ARC62371-ABflo594

Conjugate

ABflo® 594. Ex:588nm. Em:604nm.

Recommended Dilutions

FC 5 µl per 10⁶ cells in
100 µl volume

Background

This gene encodes a single-pass type I membrane glycoprotein with two Ig-like C2-type domains and an Ig-like V-type domain. This protein is one of the plasma membrane components of adherens junctions. It also serves as an entry for certain mutant strains of herpes simplex virus and pseudorabies virus, and it is involved in cell to cell spreading of these viruses. Variations in this gene have been associated with differences in the severity of multiple sclerosis. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.

Immunogen Information

Gene ID

5819

Swiss Prot

Q92692

Immunogen

Recombinant Protein corresponding to a sequence within amino acids 32-158 of human CD112/Nectin-2 (NP_001036189.1).

Synonyms

NECTIN2; CD112; HVEB; PRR2; PVRL2; PVRR2; nectin-2

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification



Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.03% proclin300,0.2% BSA,pH7.3.

Validation Data



Flow cytometry: 1×10^6 MOLT-4 cells (negative control, Left) and MCF-7 cells (Right) were surface-stained with ABflo® 594 Rabbit anti-Human CD112/Nectin-2 mAb (A24329, 5 μ l/Test, orange line) or ABflo® 594 Rabbit IgG isotype control (A23821, 5 μ l/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).

Flow cytometry: 1×10^6 MCF-7 cells were surface-stained with ABflo® 594 Rabbit IgG isotype control (A23821, 5 μ l/Test, left) or ABflo® 594 Rabbit anti-Human CD112/Nectin-2 mAb (A24329, 5 μ l/Test, right).