

HSFA2 Rabbit pAb

Catalog No.: A24390

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

Cross-Reactivity

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

817155

Swiss Prot

O80982

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ATHSFA2; heat shock transcription factor A2; T19L18.4; T19L18_4; HSFA2

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.