

PLD6 Rabbit pAb

Catalog No.: A24761

Basic Information

Observed MW

Refer to figures

Calculated MW

28kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables cardiolipin hydrolase activity and protein homodimerization activity. Involved in mitochondrial fusion. Acts upstream of or within positive regulation of mitochondrial fusion. Located in mitochondrial outer membrane.

Recommended Dilutions

WB 1:500 - 1:1000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

201164

Swiss Prot

Q8N2A8

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ZUC; PLD6

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.