

LRP4 Rabbit pAb

Catalog No.: A25225

Basic Information

Observed MW

Refer to figures

Calculated MW

212kDa

Category

Primary antibody

Applications

ELISA,IF/ICC

Cross-Reactivity

Human

Background

This gene encodes a member of the low-density lipoprotein receptor-related protein family. The encoded protein may be a regulator of Wnt signaling. Mutations in this gene are associated with Cenani-Lenz syndrome.

Recommended Dilutions

IF/ICC 1:50 - 1:200

Immunogen Information

Gene ID

4038

Swiss Prot

O75096

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1621-1720 of human LRP4 (NP_002325.2).

Synonyms

CLSS; CMS17; LRP-4; LRP10; MEGF7; SOST2; LRP4

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.

Validation Data



Immunofluorescence analysis of paraffin-embedded Human skeletal muscle tissue using LRP4 Rabbit pAb(A25225) at a dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L)(AS007) at 1:500 dilution. Blue: DAPI for nuclear staining. Perform high pressure antigen retrieval with 0.01 M citrate buffer (pH 6.0) prior to IF staining.