

# ABflo® 647 Mouse IgG 2a isotype control

Catalog No.: A25262

## Basic Information

---

**Observed MW**

**Calculated MW**

**Category**

Primary antibody

**Applications**

FC

**Cross-Reactivity**

**CloneNo number**

AMC50006

**Conjugate**

ABflo® 647. Ex:648nm. Em:664nm.

## Background

---

## Recommended Dilutions

---

**FC**                      5 µl per 10<sup>6</sup> cells in  
   100 µl volume

## Immunogen Information

---

**Gene ID**

**Swiss Prot**

**Immunogen**

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

**Synonyms**

## Contact

---

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

---

**Source**

Mouse

**Isotype**

IgG

**Purification**

Affinity purification


**Storage**

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.03% proclin300,pH7.3.

## Validation Data

---



Flow cytometry:  $1 \times 10^6$  Raji cells (negative control, left) and Human PBMC (right) were surface-stained with ABflo® 647 Mouse anti-Human CD3 mAb (A25188, 5  $\mu$ l/Test, orange line) or ABflo® 647 Mouse IgG 2a isotype control (A25262, 5  $\mu$ l/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).