

ABflo® 488 Rabbit anti-Human CD228/MFI2 mAb

Catalog No.: A25348

Basic Information

Observed MW

Calculated MW

32kDa/80kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human

CloneNo number

ARC66603-ABflo488

Conjugate

ABflo® 488. Ex:491nm. Em:516nm.

Recommended Dilutions

FC 5 µl per 10⁶ cells in
100 µl volume

Background

The protein encoded by this gene is a cell-surface glycoprotein found on melanoma cells. The protein shares sequence similarity and iron-binding properties with members of the transferrin superfamily. The importance of the iron binding function has not yet been identified. This gene resides in the same region of chromosome 3 as members of the transferrin superfamily. Alternative splicing results in two transcript variants.

Immunogen Information

Gene ID

4241

Swiss Prot

P08582

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 20-708 of human CD228/MFI2 (NP_005920.2).

Synonyms

MTf; MFI2; MTF1; CD228; MAP97

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.03% proclin300,0.2% BSA,pH7.3.

Validation Data



Flow cytometry: 1×10^6 HUVEC cells (Low Expression, left) and A-375 cells (right) were surface-stained with ABflo® 488 Rabbit anti-Human CD228/MF12 mAb (A25348, 5 μ l/Test, orange line) or ABflo® 488 Rabbit IgG isotype control (A22069, 5 μ l/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).

Flow cytometry: 1×10^6 A-375 cells were surface-stained with ABflo® 488 Rabbit IgG isotype control (A22069, 5 μ l/Test, left) or ABflo® 488 Rabbit anti-Human CD228/MF12 mAb (A25348, 5 μ l/Test, right).