

# ABflo® 594 Mouse IgG 2b isotype control

Catalog No.: A25492

## Basic Information

---

**Observed MW**

**Calculated MW**

**Category**

Primary antibody

**Applications**

FC

**Cross-Reactivity**

**CloneNo number**

AMC50008

**Conjugate**

ABflo® 594. Ex:588nm. Em:604nm.

## Background

---

## Recommended Dilutions

---

**FC**                      5 µl per 10<sup>6</sup> cells in  
   100 µl volume

## Immunogen Information

---

**Gene ID**

**Swiss Prot**

**Immunogen**

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

**Synonyms**

## Contact

---

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

---

**Source**

Mouse

**Isotype**

IgG

**Purification**

Affinity purification


**Storage**

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.03% proclin300,pH7.3.

## Validation Data

---



Flow cytometry:  $1 \times 10^6$  Karpas 299 cells were surface-stained with ABflo® 594 Mouse IgG 2b isotype control (A25492, 5  $\mu$ l/Test, left) or ABflo® 594 Mouse Anti-Human IL-23R mAb (5  $\mu$ l/Test, right).