

# PE anti-Mouse/Rat/Human FOXP3 mAb

Catalog No.: A25653

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

47kDa

### Category

Primary antibody

### Applications

FC (intra)

### Cross-Reactivity

Human, Mouse, Rat

### CloneNo number

AMC0680

### Conjugate

PE. Ex:565nm. Em:574nm.

## Recommended Dilutions

**FC (intra)** 5 µl per 10<sup>6</sup> cells in  
100 µl volume

## Background

The protein encoded by this gene is a member of the forkhead/winged-helix family of transcriptional regulators. Defects in this gene are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX), also known as X-linked autoimmunity-immunodeficiency syndrome. Alternatively spliced transcript variants encoding different isoforms have been identified.

## Immunogen Information

### Gene ID

50943

### Swiss Prot

Q9BZS1

### Immunogen

Recombinant protein

### Synonyms

JM2; AIID; IPEX; PIDX; XPID; DIETER

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Mouse

### Isotype

Mouse IgG1, κ

### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and BSA