

APC anti-Mouse TER-119 mAb

Catalog No.: A25661

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

FC

Cross-Reactivity

Mouse

CloneNo number

ATC0012

Conjugate

APC. Ex:650nm. Em:660nm.

Background

Recommended Dilutions

FC The suggested use of this reagent is $\leq 1.0 \mu\text{g}$ per 10^6 cells in 100 μl volume. It is recommended that the reagent be titrated for optimal performance for each application

Immunogen Information

Gene ID

104231

Swiss Prot

Immunogen

Day-14 fetal liver cells from a C57BL/6 mouse

Synonyms

Ter119; TER-119

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rat

Isotype

Rat IgG2b, κ

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide

Validation Data

