

# APC anti-Mouse TER-119 mAb

Catalog No.: A25661

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

### Category

Primary antibody

### Applications

FC

### Cross-Reactivity

Mouse

### CloneNo number

ATC0012

### Conjugate

APC. Ex:650nm. Em:660nm.

## Background

---

## Recommended Dilutions

---

**FC** The suggested use of this reagent is  $\leq 1.0 \mu\text{g}$  per  $10^6$  cells in 100  $\mu\text{l}$  volume. It is recommended that the reagent be titrated for optimal performance for each application

## Immunogen Information

---

### Gene ID

104231

### Swiss Prot

### Immunogen

Day-14 fetal liver cells from a C57BL/6 mouse

### Synonyms

Ter119; TER-119

## Contact

---

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

---

### Source

Rat

### Isotype

Rat IgG2b,  $\kappa$

### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide