

FITC anti-Mouse CD335 (NKp46) mAb

Catalog No.: A25672

Basic Information

Observed MW

Refer to figures

Calculated MW

37kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Mouse

CloneNo number

ATC0022

Conjugate

FITC. Ex:491nm. Em:516nm.

Background

Acts upstream of or within defense response to virus and detection of virus. Located in cell surface. Orthologous to human NCR1 (natural cytotoxicity triggering receptor 1).

Recommended Dilutions

FC The suggested use of this reagent is $\leq 1.0 \mu\text{g}$ per 10^6 cells in $100 \mu\text{l}$ volume.

Immunogen Information

Gene ID

17086

Swiss Prot

Q8C567

Immunogen

Recombinant protein


Synonyms

Ly94; NKp46; CD335

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rat

Isotype

Rat IgG2a, κ

Purification

Affinity purification

Storage

Store at 4°C . Avoid freeze / thaw cycles.

Buffer: Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide