

Marchf2 pAb

Catalog No.: A25713

Basic Information

Observed MW

Refer to figures

Calculated MW

27kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Mouse

Background

E3 ubiquitin-protein ligase that may mediate ubiquitination of TFRC and CD86, and promote their subsequent endocytosis and sorting to lysosomes via multivesicular bodies. E3 ubiquitin ligases accept ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfer the ubiquitin to targeted substrates. Together with GOPC/CAL mediates the ubiquitination and lysosomal degradation of CFTR

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

224703

Swiss Prot

Q99M02

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

MARH2

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.