

ADP-L-glycero-D-manno-heptose Rabbit mAb

Catalog No.: A26115 **Recombinant**

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

DB,ELISA

Cross-Reactivity

Independent species

CloneNo number

ARC69506

Background

Recommended Dilutions

DB 1:500 - 1:1000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Chemical compounds corresponding to ADP-L-glycero-D-manno-heptose.

Synonyms

ADP-L-glycero-beta-D-manno-heptose

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.09% Sodium azide, 0.05% BSA, 50% glycerol, pH7.3.