

Affinity Gel-conjugated Pan Phospho-Tyrosine Rabbit mAb

Catalog No.: A26203 **Recombinant**

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

IP

Cross-Reactivity

Human, Mouse, Rat, Other (Wide Range Predicted)

CloneNo number

ARC54383-AG

Conjugate

Affinity Gel

Background

Tyrosine phosphorylation (pTyr), much of which occurred on localized multiple sites, initiates cellular signaling, governs cellular functions, and its dysregulation is implicated in many diseases, especially cancers. pTyr-specific sensing is of great significance for understanding disease states and developing targeted anticancer drugs.

Recommended Dilutions

IP 0.5µg-4µg antibody for
 400µg-600µg extracts of
 whole cells

Immunogen Information

Gene ID

Swiss Prot

Immunogen


Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.09% Sodium azide, 50% glycerol, pH7.3.