

# ABflo® 647 Rabbit anti-Mouse CD94 mAb

Catalog No.: A26485

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

21kDa

### Category

Primary antibody

### Applications

FC

### Cross-Reactivity

Mouse

### CloneNo number

ARC67768

### Conjugate

ABflo® 647. Ex:648nm. Em:664nm.

## Recommended Dilutions

FC 5 µl per 10<sup>6</sup> cells in  
100 µl volume

## Background

CD94 is a 43/39 kD C-type lectin, also known as Kp43. It is present on all NK cells, NKT cells, and a subset of CD8-positive T lymphocytes in most mouse strains. CD94 is a type-II transmembrane protein with an extracellular lectin-like domain and a short cytoplasmic tail. CD94 is expressed as a disulphide-linked heterodimer with a NKG2 subunit believed to mediate signal transduction. When associated with NKG2A, the complex triggers inhibition; when associated with NKG2C, the complex triggers stimulation. The receptor complex of CD94 and NKG2 receptors bind to the ligand, Qa-1, and are thought to play a role in maintaining self-tolerance in developing NK cells.

## Immunogen Information

### Gene ID

16643

### Swiss Prot

O54707

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

CD94

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification



Affinity purification

### Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

## Validation Data



Flow cytometry: 1X10<sup>6</sup> C57BL/6 mouse splenocytes were surface-stained with ABflo® 647 Rabbit IgG isotype control (A22070,5 µl/Test,left) or ABflo® 647 Rabbit anti-Mouse CD94 mAb (A26485,5 µl/Test,right).

Flow cytometry: 1X10<sup>6</sup> C57BL/6 mouse splenocytes were surface-stained with ABflo® 488 Rabbit anti-Mouse CD161c/NK1.1 mAb (A24922,5 µl/Test) and ABflo® 647 Rabbit IgG isotype control (A22070,5 µl/Test,left) or ABflo® 647 Rabbit anti-Mouse CD94 mAb (A26485,5 µl/Test,right).