

COX6B2 Rabbit pAb

Catalog No.: A2702 **1 Publications**

Basic Information

Observed MW

11kDa

Calculated MW

11kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to be involved in oxidative phosphorylation. Located in mitochondrial crista.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

125965

Swiss Prot

Q6YFQ2

Immunogen

Recombinant protein of human COX6B2

Synonyms

CT59; COXVIB2; COX6B2

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of A-549 cells, using COX6B2 antibody (A2702).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.