

CD82 Rabbit pAb

Catalog No.: A27194

Basic Information

Observed MW

Calculated MW

30kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Mouse

Background

Predicted to be located in plasma membrane. Predicted to be integral component of plasma membrane. Is expressed in several structures, including choroid plexus; foregut; limb segment; meninges; and olfactory epithelium. Orthologous to human CD82 (CD82 molecule).

Recommended Dilutions

ELISA

Recommended starting concentration is 1 μ g/mL.

Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

12521

Swiss Prot

P40237

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

C33; IA4; Kai1; Tspan27

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.09% Sodium azide, 50% glycerol, pH7.3.