

OsRMIL1 Rabbit pAb

Catalog No.: A27288

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

WB

Cross-Reactivity

Background

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

4335680

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

RMI1; BLAP75

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.09% Sodium azide, 50% glycerol, pH7.3.